

# PROTEKT® 500 STAND

## 34500 | SIT-TO-STAND PATIENT LIFT

The Protekt® 500 lb. Sit-to-Stand patient lift has a higher weight capacity than most other competitors sit-to-stand lifts. With its 3" dual front casters it ensures safe and easy transfers. Our performance class actuator includes an emergency quick release.

- **500 lb. weight capacity.**
- **6 sling hooks allows for multiple sling options, safer patient handling and a wider range of transfers.**
- **Battery lasts twice as long as competing brands for less downtime (24 volts).**
- **Ergonomically designed foot pedal offers hands free opening of base.**
- **Dual rear locking casters.**
- **Clip free sling hookups improve safety and reduce downtime.**



- **Full parts breakdowns and modular design on lift and electronics.**
- **Advanced electronics with auxiliary up/down lift controls.**

### SPECIFICATION FOR ITEM # 34500

Weight Capacity	500 lbs.
Base Legs Length	40"
Internal Width	19"
External Width	26"
Expandable Internal Width	38"
Base Height	4.7"
Min. Lifting Height	37"
Max. Lifting Height	71"
Net Weight	119 lbs.
Shipping Weight	134 lbs.

#### WARRANTY:

Electronics    2 years  
Frame            3 years