

# Protekt® Aire 1500

Adjustable Alternating Pressure Overlay

## Protekt® Aire 1500

Adjustable Alternating Pressure Overlay



HCPCS  
E0181

### FEATURES AND BENEFITS:

- Manual air flow adjustment knob.
- Alternating Cycle Time: 10 minutes.
- Safety end flaps secures the pad to the mattress to prevent slipping.
- Lighted On/Off switch.
- Hanging bed hooks.
- Medical grade fire retardant PVC.
- Thickness: 11.5 gauge.
- Size of pad: 35.5"W x 78.5"L.
- Weight Capacity: 300 lbs.

80010

Protekt® Aire 1500 System

## Protekt® Replacement Bubble Pads



### FEATURES AND BENEFITS:

- Medical grade fire retardant PVC.
- Size of pad: 35.5"W x 78.5"L.
- Weight Capacity: 300 lbs.
- Warranty: 18 Month Non-Prorated.

80012

Bubble Pad Replacement with Flaps



270 Washington Street, Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com