

Protekt® Aire 2000

Low Air Loss/Alternating Pressure Overlay System

HCPCS
E0372



FEATURES AND BENEFITS:

PROTEKT® AIRE 2000 PUMP

- Low air loss and alternating pressure therapies in one system.
- Air flow dial to obtain custom patient comfort and pressure levels.
- Preset 10 minute alternating cycle.
- Static function for safe and easy patient transfer.
- Visual low pressure alarm.
- Adjustable hanging brackets.

80020	Protekt® Aire 2000 36"x80"x5"
80020RR	Protekt® Aire 2000 with Raised Rails

MATTRESS & COVER

- 18 laser holes keeps the patient cool and dry.
- Removable, fluid resistant, non-shear, anti-microbial nylon cover.
- Transport connector locks in air for easy patient transfers.
- CPR valve provides rapid deflation of mattress.
- Kink resistant tubing.
- Optional raised side rails available.
- Weight Capacity: 350 lbs.
- Warranty: 18 Month Non-Prorated.



270 Washington Street, Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com