

Protekt® Aire 3000

Low Air Loss/Alternating
Pressure Mattress System



HCPCS
E0277

FEATURES AND BENEFITS:

PROTEKT® AIRE 3000 PUMP

- Low air loss and alternating pressure therapies in one system.
- Air flow dial offers custom patient comfort and pressure levels.
- Static function for safe and easy patient transfer.
- Preset 10 minute alternating cycle.
- Visual low pressure alarm indicates air flow may be compromised.

MATTRESS & COVER

- 18 laser holes keeps the patient cool and dry.
- Removable, fluid resistant, non-shear, antimicrobial nylon cover.
- Pillow top feature prevents head and neck from unnecessary movement during alternating mode.
- Optional raised side rails.
- Weight Capacity: 350 lbs.
- Warranty: 18 Month Non-Prorated.

80030	Protekt® Aire 3000 36"x80"x8"
80030RR	Protekt® Aire 3000 with Raised Rails

Protekt® Aire 3500

Low Air Loss/Alternating
Pressure Mattress System with
3" Densified Fiber Support Base



HCPCS
E0277

FEATURES AND BENEFITS:

PROTEKT® AIRE 3000 PUMP

- Low air loss and alternating pressure therapies in one system.
- Air flow dial offers custom patient comfort and pressure levels.
- Static function for safe and easy patient transfer.
- Preset 10 minute alternating cycle.
- Visual low pressure alarm indicates air flow may be compromised.

MATTRESS & COVER

- 5" mattress with 3" densified fiber base provides a safety support base in case of power failure.
- 18 laser holes keeps the patient cool and dry.
- Nylon top cover is fluid resistant and quilted for added comfort.
- Pillow top feature prevents head and neck from unnecessary movement during alternating mode.
- Optional raised side rails.
- Weight Capacity: 350 lbs. • Warranty: 18 Month Non-Prorated.

83500	Protekt® Aire 3500 36"x80"x8"
83500RR	Protekt® Aire 3500 with Raised Rails

Protekt® Aire 3600

Low Air Loss/Alternating
Pressure Mattress System
with Cell-On-Cell Support Base



HCPCS
E0277

FEATURES AND BENEFITS:

PROTEKT® AIRE 3000 PUMP

- Low air loss and alternating pressure therapies in one system.
- Air flow dial offers custom patient comfort and pressure levels.
- Static function for safe and easy patient transfer.
- Preset 10 minute alternating cycle.
- Visual low pressure alarm indicates air flow may be compromised.

MATTRESS & COVER

- Unique cell-on-cell design ensures that the mattress will hold air in the lower cells in the event of a power failure.
- 18 laser holes keeps the patient cool and dry.
- Nylon top cover is fluid resistant and quilted for added comfort.
- Pillow top feature prevents head and neck from unnecessary movement during alternating mode.
- Non-skid bottom cover prevents the mattress from slipping.
- Available in 84" length. • Optional raised side rails.
- Weight Capacity: 350 lbs. • Warranty: 18 Month Non-Prorated.

83600	Protekt® Aire 3600 36"x80"x8"
83600RR	Protekt® Aire 3600 with Raised Rails