

Protekt® Aire 3600

Low Air Loss/Alternating Pressure Mattress System
with Cell-On-Cell Support Base

HCPCS
E0277



FEATURES AND BENEFITS:

PROTEKT® AIRE 3000 PUMP

- Low air loss and alternating pressure therapies in one system.
- Air flow dial offers custom patient comfort and pressure levels.
- Static function for safe and easy patient transfer.
- Preset 10 minute alternating cycle.
- Visual low pressure alarm indicates air flow may be compromised.

83600	Protekt® Aire 3600 36"x80"x8"
83600RR	Protekt® Aire 3600 with Raised Rails

MATTRESS & COVER

- Unique cell-on-cell design ensures that the mattress will hold air in the lower cells in the event of a power failure.
- 18 laser holes keeps the patient cool and dry.
- Nylon top cover is fluid resistant and quilted for added comfort.
- Pillow top feature prevents head and neck from unnecessary movement during alternating mode.
- Non-skid bottom cover prevents the mattress from slipping.
- Available in 84" length.
- Optional raised side rails.
- Weight Capacity: 350 lbs.
- Warranty: 18 Month Non-Prorated.



270 Washington Street, Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com