

Protekt® Aire 3600AB



**Low Air Loss/
Alternating
Pressure
Mattress System
with Raised Air
Bolsters and
Cell-On-Cell
Support Base**

FEATURES AND BENEFITS:

PROTEKT® AIRE 3000 PUMP

- Low air loss and alternating pressure therapies in one system.
- Air flow dial offers custom patient comfort and pressure levels.
- Static function for safe and easy patient transfer.
- Preset 10 minute alternating cycle.
- Visual low pressure alarm indicates air flow may be compromised.

MATTRESS & COVER (Available in 84" length)

- Raised side perimeter air bolsters mitigates accidental patient roll-outs, and reduces the risk of entrapment.
- Center cut-out in air bolsters enables safe ingress and egress from the mattress.
- Unique cell-on-cell design ensures that the mattress will hold air in the lower cells in the event of a power failure.
- 18 laser holes keeps the patient cool and dry.
- Stretch nylon top cover is fluid resistant and vapor permeable.
- Pillow top feature prevents head and neck from unnecessary movement during alternating mode.
- Weight Capacity: 350 lbs.
- Warranty: 18 Month Non-Prorated.

83600AB

Protekt® Aire 3600AB 36"x80"x8"

Protekt® Aire 4600DXAB



**Low Air Loss/
Alternating
Pressure
Mattress
System with
Digital Pump,
Raised Air
Bolsters and
Cell-On-Cell
Support Base**

FEATURES AND BENEFITS:

PROTEKT® AIRE 4000DX PUMP

- Extended Warranty Digital Pump with Auto Compressor Sleep Technology (ACST) shuts off compressor at desired pressure level extending pump life cycle.
- Low air loss and alternating pressure therapies in one system.
- Adjustable weight settings allows pressure to be customized according to each patient's therapeutic requirements.
- Preset 10 minute alternating cycle.
- Static function for safe and easy patient transfer.
- Audio/visual low pressure alarm indicates air flow may be compromised.

MATTRESS & COVER (Available in 84" length)

- Raised side perimeter air bolsters mitigates accidental patient roll-outs, and reduces the risk of entrapment.
- Center cut-out in air bolsters enables safe ingress and egress from the mattress.
- Unique cell-on-cell design ensures that the mattress will hold air in the lower cells in the event of a power failure.
- 18 laser holes keeps the patient cool and dry.
- Stretch nylon top cover is fluid resistant and vapor permeable.
- Static head section prevents unnecessary head/neck movement.
- Pillow top feature prevents head and neck from unnecessary movement during alternating mode.
- Weight Capacity: 400 lbs.
- Warranty: 2 Year Non-Prorated.

84600DXAB

Protekt® Aire 4600DXAB 36"x80"x8"

Protekt® Aire 6000AB



**Low Air Loss/
Alternating
Pressure Mattress
System with
Deluxe Digital
Pump, Raised
Air Bolsters and
Cell-On-Cell
Support Base**

FEATURES AND BENEFITS:

PROTEKT® AIRE 6000 PUMP

- Dual therapies; low air loss and alternating pressure.
- 4 alternating cycles 10, 15, 20 & 25 minutes allows clinicians to customize therapy.
- Seat-inflate feature helps protect the patient when they are in the fowler position.
- Auto firm mode for quick inflation.
- Static mode provides firm support for egress or ingress from mattress or side edge seating.
- 3 audible and visual alarms; low pressure, power failure and alternating failure.

MATTRESS & COVER (Available in 84" length)

- Raised side perimeter air bolsters mitigates accidental patient roll-outs, and reduces the risk of entrapment.
- Center cut-out in air bolsters enables safe ingress and egress from the mattress.
- Unique cell-on-cell design ensures that the mattress will hold air in the lower cells in the event of a power failure.
- 18 laser holes keeps the patient cool and dry.
- Stretch nylon top cover is fluid resistant and vapor permeable.
- Pillow top feature prevents head and neck from unnecessary movement during alternating mode.
- Weight Capacity: 450 lbs.
- Warranty: 2 Year Non-Prorated.

80060AB

Protekt® Aire 6000AB 36"x80"x8"