

Protekt® Aire 4000DX

Low Air Loss/Alternating
Pressure Mattress System with
Digital Pump



FEATURES AND BENEFITS:

PROTEKT® AIRE 4000DX PUMP

- Extended Warranty Digital Pump with Auto Compressor Sleep Technology (ACST) shuts off compressor at desired pressure level extending pump life cycle.
- Low air loss and alternating pressure therapies in one system.
- Adjustable weight settings allows pressure to be customized according to each patient's therapeutic requirements.
- Preset 10 minute alternating cycle.
- Static function for safe and easy patient transfer.
- Audio/visual low pressure alarm indicates airflow may be compromised.

MATTRESS & COVER

- 18 laser holes keeps the patient cool and dry.
- Nylon top cover is fluid resistant and quilted for added comfort.
- Static head section prevents unnecessary head/neck movement.
- Optional raised side rails.
- Weight Capacity: 400 lbs. • Warranty: 2 Year Non-Prorated.

| | |
|-----------|--|
| 80040DX | Protekt® Aire 4000DX 36"x80"x8" |
| 80040DXRR | Protekt® Aire 4000DX with Raised Rails |

FEATURES AND BENEFITS:

PROTEKT® AIRE 4000DX PUMP

- Extended Warranty Digital Pump with Auto Compressor Sleep Technology (ACST) shuts off compressor at desired pressure level extending pump life cycle.
- Low air loss and alternating pressure therapies in one system.
- Adjustable weight settings allows pressure to be customized according to each patient's therapeutic requirements.
- Preset 10 minute alternating cycle.
- Static function for safe and easy patient transfer.
- Audio/visual low pressure alarm indicates airflow may be compromised.

MATTRESS & COVER

- Unique cell-on-cell design ensures that the mattress will hold air in the lower cells in the event of a power failure.
- 18 laser holes keeps the patient cool and dry.
- Nylon top cover is fluid resistant and quilted for added comfort.
- Static head section prevents unnecessary head/neck movement.
- Non-skid bottom cover prevents the mattress from slipping.
- Available in 84" length. • Optional raised side rails.
- Weight Capacity: 400 lbs. • Warranty: 2 Year Non-Prorated.

| | |
|-----------|--|
| 84600DX | Protekt® Aire 4600DX 36"x80"x8" |
| 84600DXRR | Protekt® Aire 4600DX with Raised Rails |

FEATURES AND BENEFITS:

PROTEKT® AIRE 5000DX PUMP

- Extended Warranty Digital Pump with Auto Compressor Sleep Technology (ACST) shuts off compressor at desired pressure level extending pump life cycle.
- Low air loss and alternating pressure therapies in one system.
- Adjustable weight settings allows pressure to be customized according to each patient's therapeutic requirements.
- Preset 10 minute alternating cycle.
- Static function for safe and easy patient transfer.
- Audio/visual low pressure alarm indicates airflow may be compromised.

MATTRESS & COVER

- 5" mattress with 3" densified fiber base provides a safety support base in case of power failure.
- 18 laser holes keeps the patient cool and dry.
- Nylon top cover is fluid resistant and quilted for added comfort.
- Static head section prevents unnecessary head/neck movement.
- Optional raised side rails.
- Weight Capacity: 400 lbs. • Warranty: 2 Year Non-Prorated.

| | |
|-----------|--|
| 80050DX | Protekt® Aire 5000DX 36"x80"x8" |
| 80050DXRR | Protekt® Aire 5000DX with Raised Rails |

Protekt® Aire 5000DX

Low Air Loss/Alternating Pressure
Mattress System with Digital Pump
and 3" Densified Fiber Support Base

