

# Protekt® Aire 5000DX

Low Air Loss/Alternating Pressure Mattress System with Digital Pump and 3" Densified Fiber Support Base

HCPCS  
E0277



## FEATURES AND BENEFITS:

### PROTEKT® AIRE 5000DX PUMP

- Extended Warranty Digital Pump with Auto Compressor Sleep Technology (ACST) shuts off compressor at desired pressure level extending pump life cycle.
- Low air loss and alternating pressure therapies in one system.
- Adjustable weight settings allows pressure to be customized according to each patient's therapeutic requirements.
- Preset 10 minute alternating cycle.
- Static function for safe and easy patient transfer.
- Audio/visual low pressure alarm indicates airflow may be compromised.

### MATTRESS & COVER

- 5" mattress with 3" densified fiber base provides a safety support base in case of power failure.
- 18 laser holes keeps the patient cool and dry.
- Nylon top cover is fluid resistant and quilted for added comfort.
- Static head section prevents unnecessary head/neck movement.
- Optional raised side rails.
- Weight Capacity: 400 lbs.
- Warranty: 2 Year Non-Prorated.

80050DX	Protekt® Aire 5000DX 36"x80"x8"
80050DXRR	Protekt® Aire 5000DX with Raised Rails



270 Washington Street, Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com