

# Protekt® Aire 6000

Low Air Loss/Alternating Pressure Mattress System w/ Deluxe Digital Pump, Cell-On-Cell Support Base



HCPCS  
E0277

## FEATURES AND BENEFITS:

### PROTEKT® AIRE 6000 PUMP

- Dual therapies; low air loss and alternating pressure.
- 4 alternating cycles 10, 15, 20 & 25 minutes allows clinicians to customize therapy.
- Seat-inflate feature helps protect the patient when they are in the fowler position.
- Auto firm mode for quick inflation.
- Static mode provides firm support for egress or ingress from mattress or side edge seating.
- 3 audible and visual alarms; low pressure, power failure and alternating failure.

### MATTRESS & COVER

- Unique cell-on-cell design ensures that the mattress will hold air in the lower cells in the event of a power failure.
- 18 laser holes keeps the patient cool and dry.
- Nylon top cover is fluid resistant and quilted for added comfort.
- Pillow top feature prevents head and neck from unnecessary movement during alternating mode.
- Available in 84" length.
- Optional raised side rails.
- Weight Capacity: 450 lbs.
- Warranty: 2 Year Non-Prorated.

80060	Protekt® Aire 6000 36"x80"x8"
80065	Protekt® Aire 6000 with Raised Rails

# Protekt® Aire 6400

Low Air Loss/Alternating Pressure Mattress System w/ Deluxe Digital Pump



HCPCS  
E0277

## FEATURES AND BENEFITS:

### PROTEKT® AIRE 6000 PUMP

- Dual therapies; low air loss and alternating pressure.
- 4 alternating times of 10, 15, 20 & 25 minutes allows clinicians to customize therapy.
- Seat-inflate feature helps protect the patient when they are in the fowler position.
- Auto firm mode for quick inflation.
- Static mode provides firm support for egress or ingress from mattress or side edge seating.
- 3 audible and visual alarms; low pressure, power failure and alternating failure.

### MATTRESS & COVER

- 18 laser holes keeps the patient cool and dry.
- Nylon top cover is fluid resistant and quilted for added comfort.
- Pillow top feature prevents head and neck from unnecessary movement during alternating mode.
- Available in 84" length.
- Optional raised side rails.
- Weight Capacity: 450 lbs.
- Warranty: 2 Year Non-Prorated.

86400	Protekt® Aire 6400 36"x80"x8"
86400RR	Protekt® Aire 6400 with Raised Rails

# Protekt® Aire 6450

Low Air Loss/Alternating Pressure Mattress System w/ Deluxe Digital Pump and 3" Densified Fiber Support Base



HCPCS  
E0277

## FEATURES AND BENEFITS:

### PROTEKT® AIRE 6000 PUMP

- Dual therapies; low air loss and alternating pressure.
- 4 alternating times of 10, 15, 20 & 25 minutes allows clinicians to customize therapy.
- Seat-inflate feature helps protect the patient when they are in the fowler position.
- Auto firm mode for quick inflation.
- Static mode provides firm support for egress or ingress from mattress or side edge seating.
- 3 audible and visual alarms; low pressure, power failure and alternating failure.

### MATTRESS & COVER

- 5" mattress with 3" densified fiber base provides a safety support base in case of power failure.
- 18 laser holes keeps the patient cool and dry.
- Nylon top cover is fluid resistant and quilted for added comfort.
- Pillow top feature prevents head and neck from unnecessary movement during alternating mode.
- Available in 84" length.
- Optional raised side rails.
- Weight Capacity: 450 lbs.
- Warranty: 2 Year Non-Prorated.

86450	Protekt® Aire 6450 36"x80"x8"
86450RR	Protekt® Aire 6450 with Raised Rails