

Protekt® Aire 6000AB

Low Air Loss/Alternating Pressure Mattress
System with Deluxe Digital Pump, Raised Air Bolsters and
Cell-On-Cell Support Base



FEATURES AND BENEFITS:

PROTEKT® AIRE 6000 PUMP

- Dual therapies; low air loss and alternating pressure.
- 4 alternating cycles 10, 15, 20 & 25 minutes allows clinicians to customize therapy.
- Seat-inflate feature helps protect the patient when they are in the fowler position.
- Auto firm mode for quick inflation.
- Static mode provides firm support for egress or ingress from mattress or side edge seating.
- 3 audible and visual alarms; low pressure, power failure and alternating failure.

MATTRESS & COVER (Available in 84" length)

- Raised side perimeter air bolsters mitigates accidental patient roll-outs, and reduces the risk of entrapment.

- Center cut-out in air bolsters enables safe ingress and egress from the mattress.
- Unique cell-on-cell design ensures that the mattress will hold air in the lower cells in the event of a power failure.
- 18 laser holes keeps the patient cool and dry.
- Stretch nylon top cover is fluid resistant and vapor permeable.
- Pillow top feature prevents head and neck from unnecessary movement during alternating mode.
- Weight Capacity: 450 lbs.
- Warranty: 2 Year Non-Prorated.

80060AB

Protekt® Aire 6000AB 36"x80"x8"



270 Washington Street, Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com