

# Protekt® Aire 6450

Low Air Loss/Alternating Pressure Mattress System w/  
Deluxe Digital Pump and 3" Densified Fiber Support Base

HCPCS  
E0277



## FEATURES AND BENEFITS:

### PROTEKT® AIRE 6000 PUMP

- Dual therapies; low air loss and alternating pressure.
- 4 alternating times of 10, 15, 20 & 25 minutes allows clinicians to customize therapy.
- Seat-inflate feature helps protect the patient when they are in the fowler position.
- Auto firm mode for quick inflation.
- Static mode provides firm support for egress or ingress from mattress or side edge seating.
- 3 audible and visual alarms; low pressure, power failure and alternating failure.

86450	Protekt® Aire 6450 36"x80"x8"
86450RR	Protekt® Aire 6450 with Raised Rails

### MATTRESS & COVER

- 5" mattress with 3" densified fiber base provides a safety support base in case of power failure.
- 18 laser holes keeps the patient cool and dry.
- Nylon top cover is fluid resistant and quilted for added comfort.
- Pillow top feature prevents head and neck from unnecessary movement during alternating mode.
- Available in 84" length.
- Optional raised side rails.
- Weight Capacity: 450 lbs.
- Warranty: 2 Year Non-Prorated.



270 Washington Street, Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com