

Protekt® Aire Geri-Chair

Alternating Pressure Overlay System



FEATURES AND BENEFITS:

- Designed to deliver proven alternating pressure relief therapy to patients while seated in a geri-chair or recliner.
- User friendly quiet pump provides 8 liters of air per minute.
- Stretch nylon top cover is fluid resistant and vapor permeable.
- Fits standard and bariatric geriatric recliners and is easily transportable from room to room.
- Three adjustable placement straps with buckles to secure to any style geri-chair.
- Protekt® Aire 3000 pump provides alternating pressure therapy.
- Air flow dial offers custom patient comfort and pressure levels.
- Static function for safe and easy patient transfer.
- Preset 10 minute alternating cycle.
- Low pressure alarm indicates air flow may be compromised.
- Weight Capacity: 350 lbs.
- Available in 20" and 24" widths.
- Warranty: 18 Month Non-Prorated.

80018	Protekt® Aire Geri-Chair Overlay System - 20"	69"x20"x3"
80019	Protekt® Aire Geri-Chair Overlay System - 24"	69"x24"x3"
80018-M	Protekt® Aire Geri-Chair Mattress Only	69"x20"x3"
80019-M	Protekt® Aire Geri-Chair Mattress Only	69"x24"x3"
80031	Protekt® Aire Geri-Chair Pump Only	-



270 Washington Street, Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com