

Protekt® Alarms

Advanced Chair/Bed Alarm (w/nurse call)



CHAIR ALARM, BED ALARM, AND NURSE CALL!

- **Auto-Reset** - Automatically resets when pressure is reapplied on the sensor pad.
- **Visual signaling** - Top mounted flashing lights helps to verify that the monitor is armed (slow green flash) and locate the monitor when alarming (fast red flash) - perfect for nighttime use.
- **Hidden dual button reset** - Helps deter resident tampering while offering easy access for staff.
- **Floor pad compatible** - Simply switch the control on the alarm and turn it into a floor sensor monitor.
- **Nurse call port** - The ¼" plug works with most existing nurse call systems making this the ideal choice for a bed alarm monitor.
- **Meets JCAHO guidelines** - Alarm is always in the on position, making it impossible for staff to forget to turn the alarm on, and preventing the alarm from being accidentally turned off.
- **Compatibility** - Can be used with most companies bed and chair sensor pads.
- **Mounting bracket** - Easily attaches to practically all wheelchairs and beds.
- **Mounting strap** - Easily detachable strap offers the facility more mounting options.
- **Strain relief** - Avoid having to replace sensor pads because of ripped cords.
- **Volume options** - Low and high volume settings.
- **Tone options** - Different tones allow you to choose the sound that works for you.
- **Screw option** - Optional screw locks the control panel to prevent tampering.
- Rubber case protects alarm from accidental drops.
- 2 second delay option helps reduce false alarms.
- Low battery alert.
- Operates on one 9V battery (included).
- 1 Year warranty.

10260B Advanced Chair/Bed Alarm (w/nurse call)

Advanced Ultimate Alarm (all in one w/nurse call)



MAGNET, CHAIR ALARM, BED ALARM, AND NURSE CALL! EVERYTHING YOU NEED IN ONE ALARM.

- **Auto-Reset** - Automatically resets when pressure is reapplied on the sensor pad.
- **Visual Signaling** - Top mounted flashing lights helps to verify that the monitor is armed (slow green flash) and locate the monitor when alarming (fast red flash) - perfect for nighttime use.
- **Hidden dual button reset** - Helps deter resident tampering while offering easy access for staff.
- **Floor pad compatible** - Simply switch the control on the alarm and turn it into a floor sensor monitor.
- **Nurse call port** - The ¼" plug works with most existing nurse call systems making this the ideal choice for a bed alarm monitor.
- **Meets JCAHO guidelines** - Alarm is always in the on position, making it impossible for staff to forget to turn the alarm on, and preventing the alarm from being accidentally turned off.
- **Magnet option auto sensor** - When a sensor pad is plugged in the magnet is automatically deactivated, eliminating the requirement of "dummy" magnets for deactivation.
- **Compatibility** - Can be used with most companies bed and chair sensor pads.
- **Mounting bracket** - Easily attaches to practically all wheelchairs and beds.
- **Mounting strap** - Easily detachable strap offers the facility more mounting options.
- **Strain relief** - Avoid having to replace sensor pads because of ripped cords.
- **Volume options** - Low and high volume settings.
- **Tone options** - Different tones allow you to choose the sound that works for you.
- **Screw option** - Optional screw locks the control panel to prevent tampering.
- Rubber case protects alarm from accidental drops.
- 2 second delay option helps reduce false alarms.
- Low battery alert.
- Operates on one 9V battery (included).
- 1 Year warranty.

10270B Advanced Ultimate Alarm (all in one w/nurse call)