

Basic Magnet Alarm



- **Meets JCAHO guidelines** - Alarm is always in the on position, making it impossible for staff to forget to turn the alarm on, and preventing the alarm from being accidentally turned off.
- **Mounting bracket** - Easily attaches to practically all wheelchairs or beds.
- Clip attaches securely to resident's clothing.
- Reliable multi-directional magnetic pull-switch.
- Adjustable-length cord clips to resident's clothing.
- Replacement magnets available in case of misplacement.
- Operates on one 9V battery (included).
- 1 Year warranty.

10230	Basic Magnet Alarm
10230M	Basic Alarm Replacement Magnet

Basic Chair/Bed Alarm



- **Auto-Reset** - Automatically resets when pressure is reapplied on the sensor pad.
- **Volume options** - Low and high volume settings offers the caregiver the option to adjust the alarm to the required volume.
- **Meets JCAHO guidelines** - Alarm is always in the on position, making it impossible for staff to forget to turn the alarm on, and preventing the alarm from being accidentally turned off.
- **Compatibility** - Can be used with most companies bed and chair sensor pads.
- **Mounting bracket** - Easily attaches to practically all wheelchairs and beds.
- Operates on one 9V battery (included).
- 1 Year warranty.

10210	Basic Chair/Bed Alarm
-------	-----------------------

Basic Premium Magnet Alarm



The Basic Premium Magnet alarm is the most versatile magnet alarm on the market. This alarm is a standard alarm with all your basic functions that can be used in any location. Easy to set up and use, this magnet pull cord alarm is designed specifically to withstand "tugs" from any location, this patient alarm is built to last. When the patient rises, the pull cord detaches from the magnet signaling the alarm. This alarm is built with the patients safety in mind.

- **Visual signaling battery check** - Built in battery check button with LED light enables the caregiver to check battery strength at all times.
- **On/Off Switch** - Convenient on/off switch for easy deactivation.
- Magnet can be pulled from any direction.
- **Tamper-proof battery cover** - Two screws keep the battery chamber closed preventing patient tampering.
- **Versatility** - This alarm comes standard with a mounting holster, mounting strap, and extended length pull cord, allowing it to be used in a multitude of applications.
- **Mounting bracket** - Easily attaches to practically all wheelchairs and beds.
- **Mounting strap** - Easily detachable strap offers the facility more mounting options.
- **Volume options** - Low (85 dB) and high (95 dB) volume option offers the caregiver the option to adjust the alarm to the required volume.
- Cord length adjusts from 18" up to 58".
- 1 Year warranty.

10231	Basic Premium Magnet Alarm
10231M	Basic Premium Alarm Replacement Magnet