

Protekt® Anti-Abduction Knee Spreader



FEATURES AND BENEFITS:

- Comfortably holds hips in abduction while distributing pressure over thighs.
- Available with either a polyester cover or in an easy to clean vinyl.
- Strap is 46" long.
- Positions thighs 4" apart.
- Latex free.

60130	Anti-Abduction Knee Spreader - Polyester
60140	Anti-Abduction Knee Spreader - Vinyl



270 Washington Street, Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com