





Protekt® Gel Overlay



HCPCS
E0185



FEATURES AND BENEFITS:

- Optimal overlay surface for the prevention and treatment of pressure ulcers.
- Convoluted high density polyurethane foam top layer.
- 15 horizontal gel bladders prevents migration and bottoming out.
- Gel bladders contain an aqueous gel composition that provides deep immersion.
- Base foundation is a high density/high resilient foam construction which prevents bottoming out.
- Nylon top cover is fluid resistant, vapor permeable and anti-microbial.
- Bottom vinyl cover has 4 elastic corner straps that secures overlay to the mattress.

- Gel maintains room temperature at all times.
- Size: 34"x76"x3.5".
- Also available in 42", 48" and 54" widths.
- Warranty: 1 Year Non-Prorated.

90001	Gel Overlay	34"x76"x3.5"	300 lb. wt. cap.
90002	Gel Overlay Bariatric	42"x76"x4"	500 lb. wt. cap.
90003	Gel Overlay Bariatric	48"x76"x4"	550 lb. wt. cap.
90004	Gel Overlay Bariatric	54"x76"x4"	650 lb. wt. cap.



270 Washington Street, Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com