

# PROTEKT® MOVE

## 32450 | STANDING TRANSFER

The Proakt® Move is your best choice for fall prevention, reducing caregiver lifting injuries and promoting patient mobility. The Proakt® Move offers ample weight capacity at 450 lbs. and includes an adjustable base which provides clearance around common obstacles such as recliners allowing access to patients.

- **450 lb. weight capacity.**
- **Dual shin pads for extra support.**
- **Support strap available - attachment hooks for optional side support and back strap.**
- **Simple solution to compliment your electric stand assist lifters.**
- **First choice for any fall prevention program.**

### SPECIFICATION FOR ITEM # 32450

Weight Capacity	450 lbs.
Minimum Base Width	26.8"
Open Base Width	37"
Overall Height	43.3"
Minimum Seat Height	26.7"
Seat Width	20"
Overall Length	35.4"
Shin Pads	yes
Optional Support Strap	yes
Optional Seat Locks	yes

**WARRANTY:**  
Frame 3 years

