

Protekt® ProTemp Infrared Non-Contact Thermometer



The FDA approved Protekt® Pro-Temp Non-Contact Infrared Thermometer has been developed using the latest infrared technology, allowing the temporal artery (TA) temperature to be taken without touching the patient, at a distance of 2 inches away from the forehead. Precise, instantaneous, and without the need of any physical contact, the Pro-Temp is the most suitable thermometer for no risk temperature measurement. Simply point the thermometer at the patient's forehead, click "Scan" and their temperature will be instantly displayed.

FEATURES AND BENEFITS:

- Designed to take the body temperature of a person without the interference from the heat of the surrounding environment.
- Never wake the patient to take their temperature.
- Audio alarm when the temperature is too high. [$\geq 100.4^{\circ}\text{F}$ (38°C)]
- Non-touch reduces cross contamination.
- Memory for 32 previous measurements.
- 3 Color back-lit LCD screen. (Green color: Normal temperature; Orange color: Low fever; Red color: High fever).
- Three "talking" language options; English, Spanish & French.
- Data can be displayed in either Fahrenheit or Celsius.
- Body-mode temperature range: 89.6°F - 109.4°F (32°C - 42.9°C)
- Measuring distance: 1.2 in-2 in (3cm-5cm)
- Accuracy: $\pm 0.6^{\circ}\text{F}$ (0.3°C)
- Dimensions: Individual unit 6.9"(l)x1.9"(w)x2.1"(h)/ Weight: 0.2 lbs.
- Automatic shut off.
- Fast and accurate, 1 second response time.
- Hygienic, practical, and easy to use.
- No probe covers reduces costs.
- Convenient storage case included.
- Works with 2 AA batteries (included).
- FDA approved.
- Warranty: 6 Month limited.

40010

Infrared Non-Contact Thermometer



270 Washington Street, Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com