

Protekt® Safety Foam Fall Mat



FEATURES AND BENEFITS:

Manufactured from specially cured high-density impact resistant foam. The Protekt® Safety Foam Fall Mat helps reduce injuries by providing full length coverage of the bedside area.

- Bi-fold design makes storage easy, convenient and fast.
- Durable easy to clean nylon top cover is antimicrobial/antibacterial and specially designed for long term care facilities.
- The durable vinyl bottom cover mitigates the mat from slipping.
- The Protekt® Safety Foam Fall Mat is the first step in ensuring your patients safety.
- Fire retardant foam meets CA #117.
- 4 different sizes to choose from.
- Custom sizes available upon request.

51003	Protekt® Safety Foam Fall Mat - Bi-Fold	24"x66"x2"
51004	Protekt® Safety Foam Fall Mat - Bi-Fold	36"x66"x2"
51007	Protekt® Safety Foam Fall Mat - Bi-Fold	24"x66"x3"
51008	Protekt® Safety Foam Fall Mat - Bi-Fold	36"x66"x3"
51003-TF	Protekt® Safety Foam Fall Mat - Tri-Fold	24"x66"x2"
51004-TF	Protekt® Safety Foam Fall Mat - Tri-Fold	36"x66"x2"
51007-TF	Protekt® Safety Foam Fall Mat - Tri-Fold	24"x66"x3"
51008-TF	Protekt® Safety Foam Fall Mat - Tri-Fold	36"x66"x3"



270 Washington Street, Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com