

Protekt® Seat Relief

Alternating Pressure Seat Cushion System



FEATURES AND BENEFITS:

- Compact pump is rich with customizable features, including; alternating cycle times of 10, 15 or 20 min, (5) comfort settings, and static function.
- User friendly pump with rechargeable battery that can last for 12 hours with continuous operation.
- VGPT™ foam top layer offers the latest in gel-infused memory foam technology for added pressure redistribution and comfort.
- Low Battery & Low Pressure indicators ensure that the system is always working properly.
- Convenient carrying case that makes the system easily portable.
- Available in 16", 18", 20" and 22" widths.
- Weight Capacity: 300 lbs.
- Warranty: 1 Year Non-Prorated.

80120	Protekt® Seat Relief System	16"x16"x3"
80121	Protekt® Seat Relief System	18"x16"x3"
80124	Protekt® Seat Relief System	18"x18"x3"
80122	Protekt® Seat Relief System	20"x16"x3"
80126	Protekt® Seat Relief System	20"x18"x3"
80127	Protekt® Seat Relief System	22"x18"x3"



270 Washington Street, Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com