


# Protekt® Seating & Positioning



## Protekt® Foam Cushion HCPCS E2601

- Helps in the prevention of pressure sores.
- High density polyurethane foam.
- Breathable and comfortable stretch cover.
- Safety buckles keeps the cushion in place.
- 1 year warranty. • Weight Capacity: 250 lbs.

70001	16"x16"x2"	71003	18"x18"x3"
70002	18"x16"x2"	71004	20"x16"x3"
70003	20"x16"x2"	71008	20"x18"x3"
71014	14"x14"x3"	72001	16"x16"x4"
71001	16"x16"x3"	72002	18"x16"x4"
71002	18"x16"x3"	72003	20"x16"x4"



## Protekt® Ultra Cushion HCPCS E2607

- Molded high density pressure sensitive foam.
- Gel Bladder designed to eliminate migration.
- Gel infused visco foam top layer.
- Soft foam allows immersion of pelvis for comfort and stability.
- Medial and lateral side supports facilitates positioning.
- Trochanter cut outs facilitates proper pressure redistribution.
- Low shear and breathable stretch nylon cover.
- Safety buckles keeps the cushion in place.
- 18 month warranty. • Weight Capacity: 275 lbs.


77001	16"x16"x3"	77003	20"x16"x3"
77002	18"x16"x3"	77009	20"x18"x3"
77008	18"x18"x3"		



## Protekt® Supreme Bariatric Cushion HCPCS E2606

- Ultra high density pressure sensitive molded foam.
- Leg troughs promote thigh alignment and postural symmetry.
- Raised front ridge controls sliding.
- Low shear and breathable stretch nylon cover.
- Safety buckles keeps the cushion in place.
- 18 month warranty. • Weight Capacity: 500 lbs.


76005SP1	22"x16"x3"	76006SP	24"x18"x3"
76005SP	22"x18"x3"	76007SP	26"x18"x3"



## Protekt® Gel Cushion HCPCS E2603

- Helps in prevention and treatment of pressure sores.
- Gel bladder designed to eliminate migration.
- High density polyurethane foam.
- Breathable and comfortable stretch cover.
- Safety buckles keeps the cushion in place.
- 18 month warranty. • Weight Capacity: 275 lbs.

73001	16"x16"x2"	74003	18"x18"x3"
73002	18"x16"x2"	74004	20"x16"x3"
73003	20"x16"x2"	74014	20"x18"x3"
74018	14"x14"x3"	75001	16"x16"x4"
74001	16"x16"x3"	75002	18"x16"x4"
74002	18"x16"x3"	75003	20"x16"x4"



## Protekt® Foam Bariatric Cushion HCPCS E2602

- Ultra high density foam.
- Helps in the prevention of pressure sores.
- Breathable and comfortable stretch cover.
- Safety buckles keeps the cushion in place.
- 1 year warranty. • Weight Capacity: 400 lbs.

71005	22"x18"x3"	72004	22"x18"x4"
71006	24"x18"x3"	72005	24"x18"x4"
71012	24"x20"x3"	72106	24"x20"x4"
71007	26"x18"x3"	72006	26"x18"x4"
71017	26"x20"x3"	72107	26"x20"x4"
71010	28"x20"x3"	72007	28"x20"x4"
71016	30"x20"x3"	72046	30"x20"x4"



## Protekt® Ultra Bariatric Cushion HCPCS E2608

- Ultra high density pressure sensitive foam.
- Gel Bladder designed to eliminate migration.
- Gel infused visco foam top layer.
- Soft foam allows immersion of pelvis for comfort and stability.
- Medial and lateral side supports facilitates positioning.
- Trochanter cut outs facilitates proper pressure redistribution.
- Low shear and breathable stretch nylon cover.
- Safety buckles keeps the cushion in place.
- 18 month warranty. • Weight Capacity: 500 lbs.


77004	22"x18"x3"	77006	26"x18"x3"
77010	22"x20"x3"	77012	26"x20"x3"
77005	24"x18"x3"	77014	28"x20"x3"
77011	24"x20"x3"	77013	30"x20"x3"



## Protekt® Supreme Cushion HCPCS E2605

- Molded high density pressure sensitive foam.
- Leg troughs promote thigh alignment and postural symmetry.
- Raised front ridge controls sliding.
- Low shear stretch nylon cover.
- Safety buckles keeps the cushion in place.
- 18 month warranty. • Weight Capacity: 275 lbs.

76017SP	14"x14"x3"	76003SP	18"x18"x3"
76001SP	16"x16"x3"	76004SP	20"x16"x3"
76002SP	18"x16"x3"	76008SP	20"x18"x3"



## Protekt® Gel Bariatric Cushion HCPCS E2604

- Ultra high density top foam layer.
- Gel bladder designed to eliminate migration.
- Breathable and comfortable stretch cover.
- Safety buckles keeps the cushion in place.
- 18 month warranty. • Weight Capacity: 500 lbs.

74005	22"x18"x3"	75004	22"x18"x4"
74006	24"x18"x3"	75005	24"x18"x4"
74008	24"x20"x3"	75009	24"x20"x4"
74007	26"x18"x3"	75006	26"x18"x4"
74020	26"x20"x3"	75010	26"x20"x4"
74021	28"x20"x3"	75016	28"x20"x4"
74017	30"x20"x3"	75017	30"x20"x4"



## Protekt® Back Cushion HCPCS E2611/HCPCS E2612

- Improves seating posture while providing back support and comfort.
- Conforms to provide support to lumbar region.
- Stabilizing board insert prevents hammocking in sling back wheelchair.
- Low shear and breathable stretch nylon cover.
- Adjustable safety straps for added security.
- 18 month warranty.

79001	16"x17"	79005	24"x19"
79002	18"x17"	79006	26"x19"
79003	20"x19"	79009	28"x19"
79004	22"x19"	79010	30"x19"