


Protekt® Sensor Pads

Weight-sensitive pressure pads are for bed and chair monitor use. When connected, the pressure pad triggers the monitor when weight is removed from the pressure pad alerting the caregiver to respond. The alarm automatically resets when the patient returns to the safe position, reducing false alarms.


- Pads work with all major manufacturer's alarm monitors.
- All pads are constructed with durable vinyl cover, moisture resistant and easy to clean.
- Sensor pads are designed without holes to increase their water resistance and lifespan.
- Classic chair sensor pads can be used on a variety of chairs and wheelchairs.
- Slim design provides patients with a comfortable seating environment.
- Excellent sensor function warns caregivers of patients attempting to leave chair or bed earlier than most other pads.
- Large pad size helps reduce false alarms.

Classic Chair Sensor Pad	45 day	10110	7"x15"
	6 month	10111	7"x15"
	1 year	10112	7"x15"
Classic Bed Sensor Pad	45 day	10120	10"x30"
	6 month	10121	10"x30"
	1 year	10122	10"x30"
Extra Wide Bed Sensor Pad	45 day	10130	20"x30"
	6 month	10131	20"x30"
	1 year	10132	20"x30"

Chair Sensor Pad



Bed Sensor Pad



Floor Sensor Pad

- Alarm sounds when residents exit the bed and step on the pressure-sensitive pad.
- Large, heavy-duty, nonslip grips ensures that the sensor pad does not move.
- Thin edges let wheelchairs roll onto the pad, and prevent residents from accidentally tripping on the pad.
- Extra large pads measure 20" x 30" ensuring adequate coverage.

Floor Sensor Pad - Basic	1 year	10142	20"x30"
Floor Sensor Pad - Advanced	1 year	10145	20"x30"



Alarm Door Banner

- High visibility combined with a magnetically activated alarm creates an effective deterrent to wandering.
- Alarm sounds when the door is opened.
- Sold with our Basic Premium Magnet Alarm.
- Attaches to doorpost with hook and loop.
- Fits doors from 30"- 42" (W).

10161	Protekt® Alarm Door Banner System
10161B	Protekt® Door Banner Only 42"
10231	Protekt® Basic Premium Alarm