


Sit to Stand/Transfer Slings

ITEM #	DESCRIPTION
30129 Medium	Sit to stand sling - solid: medium 440 lb. weight capacity
30130 Large	Sit to stand sling - solid: large 440 lb. weight capacity
30131 X-Large	Sit to stand sling - solid: X-large 440 lb. weight capacity
30139 XX-Large	Sit to stand sling - solid: XX-large 600 lb. weight capacity
30134 Large	Transfer sling - solid: large 350 lb. weight capacity
30136 X-Large	Transfer sling - solid: X-large 440 lb. weight capacity
30669B	Buttocks Support Strap 600 lb. weight capacity
31230	Knee Strap Brace



Stand Assist Sling



Transfer Sling



Knee Strap Brace



Buttocks Support Strap

*All slings have been independently weight capacity certified by Intertek Testing Services according to ISO10535 standards.



270 Washington Street, Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com