

Protekt® Aire 6500

Low Air Loss/Alternating Pressure Mattress System
660 LBS. Weight Capacity

HCPCS
E0277



FEATURES AND BENEFITS:

- Dual compressor with 18 LPM (liter per minute).
- Alternating pressure cycles of 10, 15, 20 & 25 allows clinicians to customize pressure redistribution and skin surface interface.
- Cell-on-cell mattress design prevents patient from bottoming out in case of power failure.
- 10" air cells are constructed of Nylon/PU for added comfort and support.
- Low air loss function controls and manages moisture to prevent skin maceration.
- 3 audio/visual alarms; low pressure, power failure, and alternating failure alerts clinicians when preset pressure levels decrease.
- Seat inflate increases air flow by 20% to middle bladders to prevent bottoming out.
- Static mode allows for easy transfer in and out of bed.
- Pillow top feature prevents head and neck from unnecessary movement during alternating mode.
- Lock out feature prevents unauthorized changes to the pump settings.
- Nylon quilted top cover is fluid resistant, breathable, vapor permeable and fire retardant.
- CPR pull tab quickly releases air from mattress.
- Weight Capacity: 660 lbs.
- Warranty: 2 Year Non-Prorated.

86500	Protekt® Aire 6500 36"x80"x10"
86500RR	Protekt® Aire 6500 with Raised Rails



270 Washington Street, Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com