


# Protekt® Fiber Homecare Mattress with Cool Comfort Cover™



HCPCS  
E0272

## FEATURES AND BENEFITS:

- 100% virgin densified polyester fiber core is designed for use in Home Care.
- Protekt® Cool Comfort™ cover is constructed of polypropylene, a breakthrough healthcare textile.
- The cover is fluid resistant, breathable, 100% recyclable, odor, stain and mold resistant.
- Lighter weight than traditional innerspring mattresses.
- Lower shipping costs than innerspring mattresses.
- Meets Federal Fire Code 16 CFR part 1633.
- Weight Capacity: 300 lbs.
- Warranty: 1 Year Non-Prorated.
- Also available with a standard nylon cover.

98008-PP	Protekt® Fiber Mattress W/PP cover	35"x76"x6"
98001-PP	Protekt® Fiber Mattress W/PP cover	35"x80"x6"
98007-PP	Protekt® Fiber Mattress W/PP cover	35"x84"x6"
98008	Protekt® Fiber Mattress W/Nylon cover	35"x76"x6"
98001	Protekt® Fiber Mattress W/Nylon cover	35"x80"x6"
98007	Protekt® Fiber Mattress W/Nylon cover	35"x84"x6"
98003	Protekt® Fiber Mattress W/Nylon cover	42"x80"x6"



270 Washington Street, Mount Vernon, NY 10553

**Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005**

[www.proactivemedical.com](http://www.proactivemedical.com) • [sales@proactivemedical.com](mailto:sales@proactivemedical.com)