

Protekt® Gap Fillers Entrapment Protection

ENSURE THAT ALL GAPS AROUND THE MATTRESS ARE FILLED

Protekt® Gap Fillers are custom-made to fit any gap and protect the resident from bed entrapment. They are ideal for mitigating the danger of entrapment in Zones 3 and 7. By closing the gaps around the perimeter of the mattress, the risk of entrapment of the head, neck, chest or limbs that may occur at the headboard, footboard and rails is eliminated.

General features of the Protekt® Gap Fillers:

- Fills the gaps between the mattress and headboard/footboard/side rails.
- Help prevent head, neck, chest and limb entrapment in Zones 3 and 7.
- Easy-to-clean vinyl cover.
- Custom-made to any size ensuring a perfect fit every time.
- Meets fire rating CA #117.
- Warranty: 1 Year Non-Prorated.



Protekt® Rail Guard

- Fills the gap between the mattress and the rail.
- Attaches securely to the rail with buckles.
- Prevents entrapment along the side of the mattress.
- Works with most rails when raised or lowered.


58300	Protekt® Rail Guard	Half - 36"x6"x3"
58301	Protekt® Rail Guard	Full - 70"x6"x3"

Protekt® Mattress Extender

- Prevents entrapment along the head/footboard.
- Fills the gap between the mattress and the footboard.
- Safely accommodates taller patients.
- Optional buckles secure the extender to the head/footboard.
- Sizes in 42", 48", 54" & 60."



EXT3663	Protekt® Mattress Extender	36"x6"x3"
EXT3664	Protekt® Mattress Extender	36"x6"x4"
EXT3665	Protekt® Mattress Extender	36"x6"x5"
EXT3666	Protekt® Mattress Extender	36"x6"x6"



Protekt® Wedge Rail Guards

- Prevents entrapment along the side of the mattress.
- Fills the gap between the mattress and the rail.
- Doubles as a side rail gap filler and a patient positioning device.
- Attaches securely to the rail with buckles.

58400	Protekt® Wedge Rail Guards	Half - 36"x12" (each)
58401	Protekt® Wedge Rail Guards	Full - 70"x12" (each)

Protekt® Wall Guard

- High-density polyurethane core helps protect walls from damage caused by beds.
- Nylon cover is antimicrobial and highly fluid resistant.
- Multiple tool-free methods allow for easy attachment to most healthcare beds.
- Sizes 36", 42", 48", 54" & 60".



58601	Protekt® Wall Guard	36"x6"x3"
58602	Protekt® Wall Guard	42"x6"x3"
58603	Protekt® Wall Guard	48"x6"x3"
58604	Protekt® Wall Guard	54"x6"x3"
58605	Protekt® Wall Guard	60"x6"x3"



270 Washington Street, Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622)

Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com