


Protekt® Geri-Chair Overlay

Dual Chamber Gel-Bladder is designed to eliminate gel migration and regulate body temperature (optional).

VGPT™ foam reduces occurrence of pressure sores by channeling heat away from the body significantly faster than traditional foam, virtually eliminating perspiration.



(Pictured with gel bladder option)

FEATURES AND BENEFITS:

The Protekt® Geri-Chair Overlay is an ultra comfortable choice to reduce interface pressure while seated in a geri-chair or recliner.

- Top layer of foam features Visco Gel Polymer Technology (VGPT™) which provides rapid heat dissipation and excellent pressure redistribution.
- Seat Area - Protekt® Gel Cushion contains a center seam gel bladder which eliminates migration.
- Top Cover - Stretch fluid resistant nylon.
- Bottom Cover - Non-skid vinyl with 3 safety straps.
- Weight Capacity: 300 lbs. (19")
500 lbs. (24")
- Warranty: 1 Year Non-Prorated.

90011	Protekt® Geri-Chair Foam Overlay w/Gel	Geri-Chair Overlay w/Gel Infused Visco Memory Foam and Gel Bladder	19"x72"x3"
90011F	Protekt® Geri-Chair Foam Only Overlay	Geri-Chair Foam Overlay w/Gel Infused Visco Memory Foam	19"x72"x3"
90012	Protekt® Geri-Chair Foam Overlay w/Gel	Geri-Chair Overlay w/Gel Infused Visco Memory Foam and Gel Bladder	24"x72"x3"
90012F	Protekt® Geri-Chair Foam Only Overlay	Geri-Chair Foam Overlay w/Gel Infused Visco Memory Foam	24"x72"x3"



270 Washington Street, Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com