


Protekt® Heel Loft Pillow With VGPT™




The Protekt® Heel Loft Pillow completely off-loads pressure on the heels by comfortably elevating the feet and suspending the heels. In addition to inhibiting the development of pressure ulcers, it also alleviates lower back pain by relaxing the spine. The Visco Gel Polymer Technology (VGPT™) offers a gel-infused memory foam top layer that provides unprecedented pressure redistribution, rapid heat dissipation, and vast increases in breathability compared to traditional foam redistribution. The high-performance, breathable stretch nylon cover, is fluid resistant and antimicrobial.

- Completely removes pressure from the heels.
- Prevents the formation of pressure ulcers.
- Helps heal existing pressure ulcers.
- Optional gel bladder adds superior comfort and cooling.
- Non-skid bottom prevents sliding.
- Warranty: 1 Year Non-Prorated.

60150	Heel Loft Pillow with VGPT™	24"x16"x7"x5.5"
60150G	Heel Loft Pillow with VGPT™ with Gel Bladder	24"x16"x7"x5.5"

Protekt® Positioning Wedge



The Protekt® Positioning Wedge is used to comfortably reposition immobile patients. By reducing pressure on vulnerable areas, it assists with the prevention and management of pressure ulcers. The wedge is designed with a 30-degree lateral incline to help with proper positioning and to relieve pressure on the hip bone and coccyx. The high-performance breathable stretch nylon cover, designed for pressure relief and comfort, is fluid resistant, anti-microbial and soft to the touch.

- Recommended incline angle of 30° helps relieve pressure on hip bone and coccyx.
- High-density polyurethane foam construction.
- Breathable stretch nylon cover reduces shearing.
- Non-skid bottom prevents sliding.
- Warranty: 1 Year Non-Prorated.

50030	30° Positioning Wedge	7"x12"x18"
50031	30° Positioning Wedge	7"x12"x24"
50032	30° Positioning Wedge	7"x12"x34"

Protekt® Footrest Extender/Legrest Pad

The Protekt® Footrest Extender provides a base of support and elevates the legs, improving seat alignment, posture, and comfort while alleviating foot drop symptoms.

- Protects lower leg from contact with wheelchair.
- Prevents feet from slipping behind footrests.
- Mitigates foot drop with removable extender that raises the footrest by 1, 2 or 3 inches.
- Secures easily with quick-release straps.
- Available in sizes covering wheelchair widths of 16" – 24".
- Warranty: 1 Year Non-Prorated.



47001	Footrest Extender/Leg Rest Pad w/1" Pad	16"-20"
47002	Footrest Extender/Leg Rest Pad w/2" Pad	16"-20"
47003	Footrest Extender/Leg Rest Pad w/3" Pad	16"-20"
47004	Footrest Extender/Leg Rest Pad w/1" Pad	22"-24"
47005	Footrest Extender/Leg Rest Pad w/2" Pad	22"-24"
47006	Footrest Extender/Leg Rest Pad w/3" Pad	22"-24"

Protekt® Legrest Pad

The Protekt® Legrest Pad provides a cushion support behind the calves and prevents the feet from getting caught under the wheelchair by blocking them from slipping behind the footrest.

- Protects lower leg from contact with wheelchair.
- Prevents feet from slipping behind footrests.
- Secures easily with quick-release straps.
- Sizes covering wheelchair widths of 16" – 24".
- Warranty: 1 Year Non-Prorated.



47020	Leg Rest Pad Only	16"-20"
47024	Leg Rest Pad Only	22"-24"