

Proactive Medical Seating & Positioning Selection Guide

| Protekt Cushion | Foam | Gel | Supreme | Ultra | O2 | Foam Coccyx | Gel Coccyx | Foam Wedge | Gel Wedge | Foam Wedge w/ Pommel | Gel Wedge w/Pommel | Anti-Thrust | Seat Relief | Seat & Back Combo | Back Cushion | Adjustable Tension Back | Geri-Chair Overlay | |
|------------------------|---------------|---------------|----------------|----------------|-----------------|---------------|---------------|---------------|---------------|----------------------|--------------------|---------------|---------------|-------------------|----------------------------|-------------------------|--------------------|--------------|
| Polyurethane Foam | X | X | | | | X | X | X | X | X | X | X | | X | X | X | X | |
| Gel | | Dual Chamber | | Quad Chamber | | | Dual Chamber | | Dual Chamber | | Dual Chamber | Dual Chamber | | | | | Optional (Dual) | |
| Molded Foam | | | X | X | | | | | | | | | | | | | | |
| Gel-Infused Visco Foam | Optional | Optional | Optional | X | | | | | | | | | Insert | | | | X | |
| Air Flotation | | | | | X | | | | | | | | | | | | | |
| Alternating Pressure | | | | | | | | | | | | | X | | | | | |
| Solid Insert | | | | | | | | | | | | | | | X | | | |
| Safety Buckle | X | X | X | X | | X | X | X | X | X | X | X | X | X | X | X | X | |
| Top Cover Material | Stretch Nylon | Stretch Nylon | Stretch Nylon | Stretch Nylon | Nylon Lycra | Stretch Nylon | Stretch Nylon | Stretch Nylon | Stretch Nylon | Stretch Nylon | Stretch Nylon | Stretch Nylon | Stretch Nylon | Stretch Nylon | Stretch Nylon | Nylon | Stretch Nylon | |
| Bottom Cover Material | Vinyl | Vinyl | Non-Skid Vinyl | Non-Skid Vinyl | Non-Skid Rubber | Vinyl | Vinyl | Vinyl | Vinyl | Vinyl | Vinyl | Vinyl | Vinyl | Vinyl | Vinyl | Nylon | Non-Skid Rubber | |
| Waterproof | X | X | X | X | ~ | X | X | X | X | X | X | X | X | X | X | X | X | |
| Standard Sizes | 14"-20" | 14"-20" | 14"-20" | 16"-20" | 16"-20" | 16"-20" | 16"-20" | 16"-20" | 16"-20" | 16"-20" | 16"-20" | 16"-20" | 16"-20" | 16"-20" | 16"-20"/32"L | 16"-20" | 16"-20" | 19"-20"/72"L |
| Bariatric Sizes | 22"-32" | 22"-32" | 22"-26" | 22"-30" | 22" | 22"-30" | 22"-30" | 22"-26" | 22"-26" | 22"-26" | 22"-26" | 22"-26" | 22"-26" | 22" | 22"-26"/32"L | 22"-30" | 22"-26" | 22"-24"/72"L |
| Heights | 2" to 4" | 2" to 4" | 2.5" | 3.5" | 2" or 4" | 2" to 4" | 2" to 4" | 3"-1.5"/4"-2" | 3"-1.5"/4"-2" | 3"-1.5"/4"-2" | 3"-1.5"/4"-2" | 5"-3" | 3" | 3" | 17" or 19" | 16" | 3" | |
| Weight Capacity (lbs.) | Up to 750 | Up to 750 | Up to 750 | Up to 750 | Up to 500 | Up to 750 | Up to 750 | Up to 750 | Up to 750 | Up to 750 | Up to 750 | Up to 750 | Up to 300 | Up to 500 | ~ | ~ | Up to 500 | |
| Warranty | 1 Year | 18 Month | 18 Month | 18 Month | 5 Year | 1 Year | 1 Year | 1 Year | 1 Year | 1 Year | 1 Year | 1 Year | 1 Year | 1 Year | 18 Month | 1 Year | 1 Year | |
| HCPCS Code | E2601/E2602 | E2603/E2604 | E2605/E2606 | E2607/E2608 | E2624/E2625 | ~ | ~ | ~ | ~ | ~ | ~ | ~ | ~ | ~ | E2601/E2611 E2602/E2612 | E2611/E2612 | E2611/E2612 | ~ |

*Most competitors styles available ** Custom styles and sizes welcomed

Private label program available

For more information, please call 855-237-7622

www.ProactiveMedical.com