

Protekt® Oxygen Products



PROTEKT® OXYGEN CYLINDER D/E CART

The Protekt® Oxygen Cylinder D/E Cart is lightweight, yet durable and allows for easy and convenient transport of an oxygen tank. The frame is constructed of high-grade steel and an attractive silver vein finish. The adjustable handle accommodates a wide range of heights for individual preferences. Two large composite wheels allow for easy, stable, and safe maneuvering over most terrains. The two angled legs with anti-slip tips offers security and stability when stationary.

FEATURES AND BENEFITS:

- Accommodates one D or E size cylinder.
- Weight capacity: 80 lbs.
- Height: 38.5" – 43"
- Weight: 5 lbs.
- Warranty: Limited Lifetime.

PM CART

Protekt® Oxygen Cylinder D/E Cart – (Case quantity: 6 units)

PROTEKT® OXYGEN REGULATORS

FEATURES AND BENEFITS:

- Proactive regulators feature a protective contents gauge and a built-in pressure relief valve for safety.
- 30 micron sintered bronze inlet filter for durability.
- Meets or exceeds accuracy standards for ASTM, American National, and CGA.
- Easy-to-operate single selector knob, laser etched flow orifices.
- Body material aluminum 6061-T6, anodized green.



- CGA 870/CGA 540 inlet connections.
- Inlet pressure 300-3000 PSI.
- Outlet pressure 50 PSI.
- Warranty: 3 Year Limited.

PMREG08B	Protekt® Oxygen Regulator, 0-8 LPM Full Size, CGA 870, Barb Outlet
PMREG015B	Protekt® Oxygen Regulator, 0-15 LPM Full Size, CGA 870, Barb Outlet
PMREG015BM	Protekt® Oxygen Regulator, 0-15 LPM Mini Size, CGA 870, Barb Outlet
PMREG015D	Protekt® Oxygen Regulator, 0-15 LPM Full Size, CGA 540, Diss Outlet
PMREG04B	Protekt® Oxygen Regulator, Pediatric, 0-4 LPM Full Size, CGA 870, Barb Outlet
PMREG04D	Protekt® Oxygen Regulator, Pediatric, 0-4 LPM Full Size, CGA 540, Diss Outlet



270 Washington Street, Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com